Corning Incorporated Life Sciences
Registered

ISO 9001:2008

# **Product Description**

Catalog Number:

356407

Product Description: Corning ® Biocoat ™ 96 well plate, clear, collagen I coated, with lid

## **Component Materials:**

Collagen I	-	Rat Tail Tendon
Resin	-	Clear Polystyrene
Lid	-	Clear Polystyrene

### **Bioburden:**

This lot has been tested and found negative for Bacteria and Fungi.

### Cell Attachment and Growth Characteristics:

This lot has been tested for the attribute of Cell Attachment and Growth utilizing human HT1080 cells, an attachment-dependent mammalian cell line in a serum free media.

### Collagen I:

The plates are coated with Collagen I which is a rat tail collagen, coated via a proprietary method.

## **Storage Conditions:**

For optimal performance, the unopened product should be stored in a dry place at 4-30 degrees C (Do not freeze).

### Handling:

The product is ready to use and does not require washing of the coated surface prior to use.

### **Performance Testing:**

Each manufacturing lot is sampled and tested in accordance with Standard Operating Procedures.

Visual and microscopic examination of the product.

Visual Attributes: Bacteria & Fungi Free : Packaging:

Turbiditv test Inspection for seal and barrier integrity, accurate labeling and correct product configuration.

#### Lot Number Designation:

8 Digit Lot Number: First 3 digits - Julian Date, start of manufacturing; Next 2 digits - Year of manufacture; Last 3 digits - Batch identification.

Rev No: 1