

sigma-aldrich.com

3050 Spruce Street, Saint Louis, MO 63103, USA Website: www.sigmaaldrich.com Email USA: techserv@sial.com Outside USA: eurtechserv@sial.com

## **Product Specification**

Product Name: Bovine Serum Albumin - pH <5.0, powder

Product Number:	A8531
CAS Number:	9048-46-8

Storage Temperature: -20 ℃

TEST	Specification
Appearance (Color)	White to Tan
Appearance (Color)	Powder
Solubility (color)	Colorless to Tan
Solubility (tubidity)	Clear
Solubility (Form)	Solution
Suitability	Suitable
Suitable for use as a gel filtration marker per Sigma Technica	al Bulletin GF-3 (MW approx. 66,000 Da)
VSV and BT Virus	None Detected
Miscellaneous Supplier Data	Conforms
Inactivation Process	
pH not more than 5.0 for at least 2 hours	
Chapter 4(D) - pH	
Product meets European Union requirements for treated techn	nical blood products. Form of treatment pH < /= $5.0$ for > /= $2$ hours.

Specification: PRD.2.ZQ5.10000010109

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.