



SLS Select Autoclave Bag, 840 x 635mm

Code: AUT1404



Dimension detail	
Product width (mm)	635mm ±10mm
Product length (mm)	840mm ±10mm
Lip (mm)	n/a
Lateral seal position	n/a
Handle position (mm)	n/a
Perforations	n/a

Print detail	
Print PDF & position	84137
Print description	n/a
Inks required	n/a
Plate type	n/a
Repeat (mm)	635

Post seal detail	
Release tape (mm)	n/a
Top glue line (mm)	n/a
2nd glue line (mm)	n/a
Glue gap (mm)	n/a
Glue type	n/a

Material detail	
Polymer	HDPE
Master batch	n/a
Additives	n/a
Film thickness (Micron)	35
Tolerance	±10 % average ±20 % spot

Packaging detail	
Pack quantity	100
Packing material	Carton
Label type	Standard
Pallet quantity	100
Pallet type (mm)	Standard
Transit packaging	Standard

United Kingdom:
slsinfo@scientific-labs.com

International:
slsexport@scientific-labs.com

Northern Ireland:
NI.sales@scientific-labs.com

Republic of Ireland:
ROI.sales@scientific-labs.com

