

Cole-Parmer®

Cole-Parmer® Electric Bunsen Mantle - BB-200 & CN-200 Analogue Controller

- Conical shaped heating element directs radiant
- Top cowl deflects heat away from your hand
- Air circulation from the vented housing keeps the base cool enough to hold during operation
- Burner consumes only 400W of power
- Recommended controller is the CN-200



Bunsen Burner BB-200

Analogue Controller CN-200



Bunsen Burner BB-200

Cole-Parmer® Electric Bunsen Mantle - BB-200 & CN-200 Analogue Controller

Cole-Parmer® Electric Bunsen BB-200, combines the advantages of a regular gas burner with the clean, easy operation of our Electromantles. It is corrosion-resistant and has a cool-to-the-touch base.

Radiation from the heater is directed upwards to a focal point. The Electric Bunsen is ideal for heating test tubes, crucibles, small flasks and beakers, regardless of their shape.

CN-200 Analogue Controller

The CN-200 analogue controller has been designed to provide a comprehensive answer to control the heating of resistive loads such as Electric Bunsens, Electromantles, Heating Tapes and Cords for bench top operation delivering up to a maximum of 800 Watts. A rod support clamp is provided at the rear of the controller to take a standard 12.5mm diameter rod.

An accessory extension mains lead is available where remote operation is required e.g. in a fume extraction unit.

Technical Specification

Specification	BB-200
Flask capacity	25ml
Material	Stainless-steel
Max temperature range	800°C to 1000°C
Shipping weight, kg	0.5
Dimensions (d x h), mm	120 x177

Ordering Information

Description (without controller)	Ordering Number	Series No.	Model No.	Legacy Sku.
Electric Bunsen & MC5 Controller 230V UK & EU	36135-02	BB-200	BB-200-CN	BA6101/C
Electric Bunsen & MC5 Controller 115V US	36135-03	BB-200	BB-200-CN-115	BA6101/CX1
Electric Bunsen 230V UK & EU	36135-01	BB-200	BB-200	BA6101
Electric Bunsen 115V US	36135-04	BB-200	BB-200-115	BA6101X1

coleparmer.com

USA:
+1.800.323.4340
+1.847.549.7600

Canada: +1.800.363.5900
China: +86.21.5109.9909
France: +33 (0) 1486 37800
Germany: +49 (0) 9377 92030

India: +1.800.266.1244
Italy: +39 (0) 1313 89513
UK: +44 (0) 1480 272279
All other countries: +1.847.549.7600

Cole-Parmer®
essentials

an antylia scientific company

Find out more at antylia.com