

Routine Laboratory Thermoregulators

Baths and Thermoregulators

The Techne® Tempette clip-on thermoregulators are designed for use in combination with the Techne® unheated water baths or any other suitable laboratory vessels. They will heat, circulate and safely control the temperature of the liquid in the bath within precise limits suitable for routine laboratory applications.

TE-10A Tempette

- Temperature range of -20°C* to 95°C
- Excellent temperature stability: $\pm 0.01^\circ\text{C}$ at 40°C
- Simple to use analogue control
- Suitable for most routine laboratory applications
- User adjustable over-temperature cut-out for unbeatable safety



TE-10A

TE-10D Tempette

- Temperature range of -40°C* to 120°C
- Excellent temperature stability: $\pm 0.01^\circ\text{C}$ at 40°C
- 4 digit setting with a bright LED digital temperature display
- Suitable for most routine laboratory applications
- User adjustable over-temperature cut-out
- Low liquid level cut-out as standard



TE-10D