High Powered Thermoregulators

Baths and Thermoregulators

Technical Specification

Specifications to DIN 12876	TE-20D	TE-20HT
Temperature range*	-20°C to +95°C	-40°C to +120°C
Temperature range*	-40°C to +200°C	-40°C to +250°C
Temperature selection	Digital	Digital
Temperature stability using water @ 40°C	±0.005°C	±0.005°C
Method of control	PID	PID
Temperature sensor	PRT	PRT
Adjustable over-temperature cut-out	Yes	Yes
Low liquid level cut-out	Yes	Yes
PC Interface	RS232	RS232
Heating/Pumping		
Pump capacity litres/minute	10	Internal circulation
Pump capacity (mbar)	145	
Nominal heater power at 120V (W)	1500	1500
Nominal heater power at 240V (W)	1800	1800
Cooling coil	No	Option
Extension below base, mm	145	145
Dimensions (L x W x H), mm	237 x 124 x 260	237 x 124 x 260
Shipping Weight, kg	4.0	4.0

^{*} Refrigeration or cooling coil required for below ambient cooling (see Techne Flow and Dip Coolers and the cooling coil). The TU-20HT can only be used with the Dip Coolers

Ordering Information

	Product codes			
Description	230V	120V	115V	100V
TU-20D, advanced thermoregulator with RS232, -40°C to 200°C, (supplied with clamp)	FTU20DDC	n/a	FTU20DPC	n/a
TU-20HT, advanced high temperature thermoregulator with RS232, -40°C to 250°C, (supplied with clamp)	FTU20HDC	FTU20HPC	n/a	FTU20HYC

Techne Catalogue Page 47