

# IKA

designed for scientists



## DB 4.6

/// Data Sheet

Used for round bottom tubes (20 mm)

Pore size: Ø 21.0 mm

No. of holes: 8

Hole depth: 48.4 mm

Dimensions (W x D x H): 95 x 76 x 51 mm

[www.ika.com](http://www.ika.com)

Subject to technical changes





designed for scientists

### Technical Data

Dimensions (W x H x D) [mm]	95 x 51 x 76
Weight [kg]	0.667
Permissible ambient temperature [°C]	5 - 40