

IKA

designed for scientists



DB 6.1

/// Data Sheet

Used for individual PCR tube strips (0.2 ml tubes)

Pore size: Ø 6.4 mm

No. of holes: 80

Hole depth: 15.5 mm

Dimensions (W x D x H): 95 x 76 x 51 mm



designed for scientists

Technical Data

| | |
|--------------------------------------|--------------|
| Dimensions (W x H x D) [mm] | 95 x 51 x 76 |
| Weight [kg] | 0.93 |
| Permissible ambient temperature [°C] | 5 - 40 |