

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.com
Email USA: techserv@sial.com
Outside USA: eurtechserv@sial.com

## **Product Specification**

Product Name:

OptiPrep™ Density Gradient Medium - used for cell and subcellular organelle isolation

Product Number: D1556

MDL: MFCD00867965

TEST Specification

Bulk Density Endotoxin Level

Sterility Trademark

OptiPrep is a trademark of Alere Technologies AS

Specification: PRD.1.ZQ5.10000033831

1.319 - 1.321 g/ml

≤ 1 EU/ml Pass

ass

Confirmed

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.