

1- glass door tall safety cabinet with filtration

COMPLIANCE

NF X 15-211, Electric fan in compliance with CE specifications

ADVANTAGES

- CE certified
 - NF X 15-211 compliant
 - Storage of flammables, corrosives and toxics
- BUT ALSO...**
- Perforated shelves
 - Lighting control allowing to survey the proper working of the fan
 - Front control window to see immediately the filter
 - Compliant to the Public Health Code articles R. 5132-66 and R. 5132-68
 - INRS recommendations compliant
 - Normalized labelling provided with the equipment
 - Available with different configurations : 1, 2, 3 or 4 compartments
 - Fan in compliance with CE specifications
 - Active coal filters adsorption capacity and filtration system certified NF X 15-211
 - Stop 99,97% of the particles with a diameter bigger than 0.3 µm
 - No modification of building's structure required
 - Active coal filter tested to the standard NF X 15-211
 - Easy replacement of the active coal filter



TECHNICAL FEATURES

Model	AF1C	Door closure	manual closing
Compliance	NF X 15-211, Electric fan in compliance with CE specifications	Shipping dimensions HxWxD (mm)	2100 x 1200 x 800
Cabinet materials	Melamine wood structure, Double-walled glass epoxy	Ventilation	Lateral inlet, outlet at the top of the cabinet
Doors color	white RAL 9010	Number of compartments	1
Number of doors	1	Containment sump	1
Doors type	swing doors	HS Code	115.00
Hinges	right		94032080

COMPLEMENTARY PRODUCTS

- Kit to use the fume hood without activated carbon filter - for cabinets range 12.X since 2008 (Réf : CAR300)
- Activated carbon filter for organic and corrosive vapors - for cabinets range 12.X since 2008 (Réf : CORG301)
- Activated carbon filter for formaldehyde vapors - for cabinets range 12.X since 2008 (Réf : FOR300)
- Activated carbon filter for organic vapors - for cabinets range 12.X since 2008 (Réf : ORG300)

OPTIONS

- Shelf for cabinets for range 12.X (Réf : ESAF2X)

Model	External dimensions HxWxD (mm)	Internal dimensions HxWxD (mm)	Weight (kg)
AF1C	1900 x 800 x 545	1500 x 760 x 500	115.00