



designed for scientists



R 104 Stand

/// Data Sheet

Small stand for T 10 basic and T 10 standard.

Diameter of support rod: 10 mm

Dimensions (W x D): 242 x 355 mm

Height: 370 mm

Max. load: 0,7 kg



designed for scientists

Technical Data

Diameter [mm]	10
Dynamic load [kg]	0.7
Dimensions (W x H x D) [mm]	242 x 370 x 355
Weight [kg]	0.878
Permissible ambient temperature [°C]	5 - 40