

Lonza Walkersville, Inc.
Walkersville, MD 21793-0127 USA
U.S. Scientific Support: 800 521 0390
scientific.support@lonza.com
EU/ROW Scientific Support: +32 87 321 611
scientific.support.eu@lonza.com
Document # TS-17-836-4 09/11
www.lonza.com
© 2012 Lonza Walkersville, Inc.

Amphotericin B

Product Use

Amphotericin B is intended for the control of yeasts and fungi in tissue culture.

Specifications

Sterility	pH	Osmolality (mOsm)
Negative	FIO	FIO

Toxicity	Endotoxin (EU/ml)
FIO	N/A

Storage

-10°C to -20°C

Ordering Information

Cat. No.	Description	Size
17-836R	Contains 250 ug/ml Amphotericin B	20 ml
17-836E	Contains 250 ug/ml Amphotericin B	100 ml

Product Use Statement

THESE PRODUCTS ARE FOR RESEARCH USE ONLY. Not approved for human or veterinary use, for application to humans or animals, or for use in clinical or *in vitro* procedures.

All trademarks herein are marks of Lonza Group or its subsidiaries.