

Phosphate Buffered Saline  
w/o Calcium, Magnesium, 17-517  
or Phenol Red (10X)



Description	CAS #	Chemical Formula	Concentration		Molarity	
			g/L	mg/L	M	mM
Potassium Phosphate Monobasic Anhydrous	7778-77-0	KH <sub>2</sub> PO <sub>4</sub>	1.440	1.440E+03	1.058E-02	10.580
Sodium Chloride	7647-14-5	NaCl	90.000	9.000E+04	1.540	1.540E+03
Sodium Phosphate, Dibasic Anhydrous	7558-79-4	Na <sub>2</sub> HPO <sub>4</sub>	7.950	7.950E+03	5.599E-02	55.986