




SMZ171 Series

# SMZ-171-BP





SMZ161 Series

**SMZ-171-BP**

|                                |  |
|--------------------------------|--|
| <b>Optical system</b>          | Greenough                                      |
| <b>Observation tube</b>        | Binocular head                                 |
| <b>Inclination</b>             | 45° inclined                                   |
| <b>Interpupillary distance</b> | 52-75mm  |
| <b>Diopter adjustment</b>      | On both eyepieces, +/- 5 diopter               |
| <b>Eyepieces</b>               | Widefield N-WF10X/23mm with diopter adjustment |
| <b>Objectives system</b>       | Zoom, ratio 6.7:1                              |
| <b>Objectives</b>              | 0.75X-5X                                       |
| <b>Working distance</b>        | 110mm  |
| <b>Stand type</b>              | Pole type stand without illumination           |
| <b>Head holder</b>             | For Ø32mm column and Ø76mm head                |
| <b>Focus mechanism</b>         | Coarse focusing system with tension adjustment |
| <b>Focusing stroke</b>         | 49mm   |
| <b>Accessories included</b>    | Dust cover, black/white stage plate            |
| <b>Dimensions LxWxH</b>        | 303x239x370mm                                  |
| <b>Base</b>                    | 303x239mm                                      |
| <b>Column/Arm</b>              | Height 290mm, Ø32mm                            |
| <b>Net weight</b>              | 5.45kg   |