# Spark-Stack<sup>TM</sup> – enhances throughput



### Spark-Stack plate stacker for walkaway automation of Tecan Spark® readers

Enhance the productivity of your lab with this compact, reliable and field-upgradeable plate stacker. Spark-Stack takes care of typical process steps, such as loading, unloading and re-stacking of non-lidded 6- to 1,536-well plates. In addition, it enables benchtop automation of light-sensitive assays, using plate magazines with protective covers.

#### **Applications:**

- Light-sensitive assays, such as AlphaLISA®, AlphaScreen® and AlphaPlex™
- Incubate-and-read assays, such as HTRF® and LanthaScreen®
- Batch-type cell-based assays, such as CellTiter-Glo® Cell Viability Assay
- Colorimetric, fluorometric and luminescence assays



- Dark covers protect light-sensitive assay plates inside the plate magazines from ambient lighting
- Enables benchtop automation of a range of light-sensitive assays, such as AlphaLISA, AlphaScreen, AlphaPlex, CellTiter-Glo, and many more

- Faster time to results thanks to overnight runs with software-controlled pre-incubation and measurement steps
- Remote monitoring of the stacker run using the Tecan CNS app and the SparkControl Magellan™ data analysis software
- Check the status of the run directly on your smartphone and benefit from true walkaway automation



## www.tecan.com/spark-stack

**Australia** +61 3 9647 4100 **Austria** +43 62 46 89 33 **Belgium** +32 15 42 13 19 **China** +86 21 220 63 206 **Denmark** +45 70 23 44 50 **France** +33 4 72 76 04 80 **Germany** +49 79 51 94 170 **Italy** +39 02 92 44 790 **Japan** +81 44 556 73 11 **Netherlands** +31 18 34 48 17 4 **Singapore** +65 644 41 886 **Spain** +34 93 490 01 74 **Sweden** +46 8 750 39 40 **Switzerland** +41 44 922 89 22 **UK** +44 118 9300 300 **USA** +1 919 361 5200 **Other countries** +43 62 46 89 33

#### For research use only.

© 2018 Tecan trading AG, Switzerland, all rights reserved. For disclaimer and trademarks please visit www.tecan.com

