

# IKA

designed for scientists



## MT 40.10 Disposable grinding chamber

/// Data Sheet

The disposable grinding chamber with a stainless steel beater reduces soft, medium, hard and brittle materials with a Mohs hardness of up to 5 (manganese or apatite: 5 Mohs). The chamber is made of transparent plastic so grinding tests can be observed at any time.

Useful volume: 40 ml

Quantity: 10 pieces/pack

Material grinding chamber: PP



designed for scientists

### Technical Data

Weight [kg]	0.394
Permissible ambient temperature [°C]	5 - 40