

IKA

designed for scientists



MT 40.100 Disposable grinding chamber

/// Data Sheet

The disposable grinding chamber with a stainless steel beater reduces soft, medium, hard and brittle materials with a Mohs hardness of up to 5 (manganese or apatite: 5 Mohs). The chamber is made of transparent plastic so grinding tests can be observed at any time

Useful volume: 40 ml

Quantity: 100 pieces/pack

Material grinding chamber: PP

www.ika.com

Subject to technical changes





designed for scientists

Technical Data

Weight [kg]	3.767
Permissible ambient temperature [°C]	5 - 40