



designed for scientists



MultiDrive BL 2000 Blender vessel

/// Data Sheet

BL 2000 is a vessel with mixing knife and is particularly useful in combination with IKA MultiDrive basic for wet grinding of a variety of applications, such as meat, sausage, soy beans, silage, vegetables.

All parts can be cleaned in a dishwasher. The knife is easy to dismantle for intensive cleaning.

The scope of delivery includes a standard mixing knife made of stainless steel.



designed for scientists

Technical Data

Process type	batch
Operating principle	blending
Speed range [rpm]	3000 - 20000
Circumferential speed max. [m/s]	78.5
Useable volume max. [ml]	2000
Feed grain size max. [mm]	30
Material beater/cutter	stainless steel 1.4034
Material (other)	stainless-steel, PCTG, PTFE, silicone
Power-ON time [min]	5
Power-OFF time [min]	10
FDA conform	yes
Multi use	yes
Dishwasher cleanable	yes
Beater/Knife diameter [mm]	75
Dimensions (W x H x D) [mm]	180 x 280 x 220
Weight [kg]	1.25
Permissible ambient temperature [°C]	5 - 40
Permissible relative humidity [%]	80

