

IKA

designed for scientists



MultiDrive MI 250.3 Beater, tungsten carbide

/// Data Sheet

Beater, tungsten carbide

Suitable for harder materials up to Mohs hardness 9

Compatible with IKA MultiDrive MI 250.

www.ika.com

Subject to technical changes





designed for scientists

Technical Data

Process type	batch
Operating principle	impact
Feed hardness max. [Mohs]	9
Material beater/cutter	Hard metal (WC)
FDA conform	yes
Multi use	yes
Beater/Knife diameter [mm]	74
Permissible ambient temperature [°C]	5 - 40