

# MULTIPLE STIRRERS WITH INTERNAL CONTROL

Designed for reliable mass screenings as well as chemical, biotechnical and medical applications of all kinds, developed for the hard daily lab use and uncompromising continuous operation.

## MIX 6 / MIX 15

Ultra-flat magnetic stirrer with 6 stirring positions, suitable for 6 x 1,000 ml Erlenmeyer flasks, or with 15 stirring positions, suitable for 15 x 250 ml beaker glasses (tall form).

- 🌀 Maintenance- and wear-free
- 🌀 For different vessel sizes
- 🌀 Outstanding long-life quality
- 🌀 Developed and made in Germany

*magnetic motion*



| Type                           | MIX 6                                 | MIX 15 |
|--------------------------------|---------------------------------------|--------|
| Stirring points                | 6                                     | 15     |
| Stirring point distance        | 130 mm                                | 65 mm  |
| Stirring volume/stirring point | 1 - 3,000 ml                          |        |
| Speed range                    | 100 - 2,000 rpm                       |        |
| Power setting                  | 5/10/15/20 W (4-step)                 |        |
| Material housing               | Stainless steel                       |        |
| Permitted operation conditions | -10 °C up to +50 °C (at 80% humidity) |        |
| Electrical data                | 100-240 V / 50-60 Hz / 1.5 A          |        |
| Dimensions (WxDxH)             | 245 x 425 x 38 mm                     |        |
| Protection category            | IP64                                  |        |
| Weight (gross)                 | approx. 8.8 kg                        |        |
| Order no.                      | 30006                                 | 30015  |

3 years warranty on material and manufacturing.  
We reserve the right to make technical alterations. We do not assume liability for errors in printing or typing.

6 or 15 stirring positions, 100% maintenance-free and wear-free by inductive 2mag-Magnetic-Drive concept, power reserves up to 3,000 ml stirring volume per stirring point, extremely wide speed range of 100 to 2,000 rpm, 100% synchronized speed, jerk-free stirring even at low speeds, 4-step power setting of the stirrer, high power setting for viscous liquids and reduced stirring power for long-term operation without any heating effects caused by the stirrer.

Clear digital display for settings of stirrer speed and power, SoftStart procedure for reliable catching/centering and safe acceleration of the stirring bar, QuickSet for quick setting of start- and maximum speed, Automatic memory function to save the last settings, outstanding long-life cycle of the hermetically sealed drive system, resistant stainless steel housing for easy and fast cleaning..

### Accessories

ASTEROID 40 - Stirring bar  
Stirring bar remover PTFE

