CONTROL UNITS FOR MIXdrive STIRRING DRIVES

Designed for 2mag MIXdrive stirring drives. For chemical, biotechnical and medical applications of all kinds.



MIXcontrol eco

Price-optimized control unit for small medium volumes and the connection of 1 stirring drive MIXdrive.

- ₱ Price-optimized
- № 120 up to 1,200 rpm
- **№** SoftStart





MIXcontrol eco
120 - 1,200 rpm
10 W
-
Plastic
0 °C up to +40 °C (at 80% humidity)
100-240 V / 50-60 Hz / 0.75 A
65 x 50 x 45 mm
IP20
approx. 0.5 kg
90100

Individual units and dimensions on request.

3 years warranty on material and manufacturing.

We reserve the right to make technical alterations. We do not assume liability for errors in printing or typeing.

Control unit with 10 watts stirring power, speed range from 100 up to 1,200 rpm, SoftStart for reliable catching/centering and safe accelerating of the stirring bar.

Compatible stirring drives

MIXdrive 1

MIXdrive 1 eco

MIXdrive 1 XS

cuvetteMIXdrive 1

MIXdrive 6

MIXdrive 15

MIXdrive 60

20 www.2mag.de