

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaldrich.com](http://www.sigmaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:

JumpStart™ Taq ReadyMix™ - Complete optimized reagent for hot-start PCR at 2X concentration

Product Number: **P2893**

Storage Temperature: -20 °C

TEST	Specification
DNase, Exonuclease Detection	None Detected
RNAse Detection	None Detected
JumpStart Taq Antibody	Pass
Reduces non-specific product formation.	
Suitable for Quantitative PCR	Suitable

Specification: PRD.0.ZQ5.10000063776

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.