

NOT TO SCALE



127.76 ±0.25

103.50 ±0.25

12.13
±0.25

4.50 typ

14.35
±0.25



SECTION X-X



WWW.4ti.co.uk

TITLE:

VISION PLATE 384

DATE:

22/3/07

SHEET No.:

1 OF 1

ISSUE:

A

All dimensions in millimetres

This drawing and/or any technical information relating to it is issued in strict confidence and may not be reproduced or used for manufacture without written permission. Upon receipt of which it is the users responsibility to seek confirmation of any dimensions and production tolerances.