



## Countertop Laboratory Refrigerator Solid Door 82L LR207C

LR207C

- ✓ 2°C to 10°C temperature range
- ✓ External digital controller and temperature display
- ✓ Fitted lock and 2 keys
- ✓ Spark free interior
- ✓ Audible and visual temperature alarms

## SPECIFICATION

Dimensions	502mm W x 660mm H x 540mm D
Weight	27.5kg
Product name	Countertop Laboratory Refrigerator Solid Door 82L LR207C
GDPA code	LR207C
SKU code	UK 444441799, EU 444442345, AU 444442346
External dimensions	502mm W x 660mm H x 540mm D
Internal dimensions	415mm W x 540mm H x 340mm D
Packed dimensions	555mm W x 677mm H x 580mm D
Product weight	27.5kg
Packed weight	30.5kg
Gross capacity	82L
Net capacity	67L

External finish	White
Colour	White/green
Internal finish	Plastic
Type of plastic	SRL800 (HIPS/ABS mix)
Door type	Solid
Door style	Hinged, reversible
Door open depth	976mm
Replacement door seal	N
Open door sensor	N
Lock	Y
Adjustable feet	N
Temperature range	2°C to 10°C
Temperature display on panel	Y
Temperature alarm	Y
Number of shelves	2
Shelf material	Wire
Number of shelf positions	2
Shelf size	438mm W x 6mm H x 274mm D
Interior light	N
Auto or manual defrost	Auto
Fan	N
Controller type	Digital, external
Climate class	N
Refrigerant type	R134a
CFC free	Y
HFC free	N
HCFC free	Y

UK 3 pin plug	Y
Voltage supply	230
Spark free	Y
Auxiliary output	N
Insulation thickness	35mm
Insulation type	Cyclopentane
Noise level	<43dB
Energy consumption	0.77 kW/24h