

KIMTECH™

Kimtech™ A7 Laboratory Coats



Breathable
multi-layered lab coats

Chemical and abrasion resistant
with cloth-like feel

Individually
packed

Kimtech™ A7 Laboratory Coats are disposable coats designed to maximise safety, comfort and compliance.

They are designed to protect the wearer from low risk liquid chemicals in the form of sprays, aerosols and light splashes. This contamination control is ensured by excellent particle and bacterial filtration, tear-resistant construction, and extensive testing to ensure regulatory compliance.

The coats are made of non-woven SMS, a secure and durable fabric with a cloth-like feel. Triple overlapped seams boost strength and durability and the breathable material decreases the risk of heat stress to enable long periods of uninterrupted wear.

The Kimtech™ A7 disposable laboratory coats are latex and silicone free, individually packed, and available in sizes S to 3XL.

Kimtech™ A7 Laboratory Coats



Key Features

- › Made from breathable and abrasion-resistant layers, the laboratory coats provide enhanced protection against chemical splash, dry particulates and a variety of water-based liquids
- › Single use, disposable garment
- › Features triple stitched and overlapped seams for added strength and durability

Assured Compliance

- › PPE Cat III according to Regulation (EU) 2016/425
- › Type 6 [PB] EN 13034 limited chemical splash protection certified
- › EN1149-1 Antistatic clothing¹

¹ This laboratory coat is antistatically treated and offers electrostatic protection according to EN1149-5:2018 when properly grounded.

Quality Standards

- › Manufactured in accordance with Quality System ISO 9001



EN 13034:2005
+A1:2009



EN 1149-5:
2018

CE 0598

Technical Data		
Fabric Test	Test Method	Result
Abrasion resistance	EN 530 M2	Class 2 (visual)
Trapezoidal tear resistance	ISO 9073-4	Class 2
Tensile strength	EN ISO 13934-1	Class 1
Puncture resistance	EN 863	Class 1
Resistance to ignition	EN 13274-4 M3	PASS
Seam strength	EN ISO 13935-2	Class 3

Kimtech™ A7 Laboratory Coats			
Code	Colour	Size	Total/Case
96700	White	S	15
96710	White	M	15
96720	White	L	15
96730	White	XL	15
96740	White	2XL	15
96750	White	3XL	15