

## 2-door large sheet steel safety cabinet

### COMPLIANCE

### ADVANTAGES

- INRS recommendations compliant
- Storage of flammables, corrosives and toxics
- Compliant to the Public Health Code

### BUT ALSO...

- Compliant to the Public Health Code articles R. 5132-66 and R. 5132-68
- INRS recommendations compliant
- Normalized labelling provided with the equipment
- Holes for the connection to the ventilation
- Allows the storage of your chemicals safely
- Affordable equipment to store your chemicals



### TECHNICAL FEATURES

Model	<b>AL305</b>	Ventilation	<b>Both inlet and outlet</b>
Cabinet materials	<b>8/10th epoxy painted steel structure</b>	Openings diameter (mm)	<b>100</b>
Doors color	<b>grey RAL 7035</b>	Number of compartments	<b>1</b>
Number of doors	<b>2</b>	Containment sump	<b>1</b>
Doors type	<b>swing doors</b>	Grating	<b>Option</b>
Hinges	<b>right &amp; left</b>	HS Code	<b>94032080</b>
Door closure	<b>manual closing</b>		
Shipping dimensions HxWxD (mm)	<b>2000 x 1200 x 800</b>		

### ADDITIONAL ACCESSORIES RECOMMENDED - FILTRATION OR VENTILATION BOXES

- Filtration box with activated carbon filter for organic vapors (Réf : CDF-A)
- Filtration box with activated carbon filter for organic and corrosive vapors (Réf : CDF-ACORG)
- Filtration box built in PVC with activated carbon filter for organic and corrosive vapors (Réf : CDF-P)
- Ventilation box for a connection to the outside (Réf : CDV-A)
- Recirculating air box Labopur - NF X 15-211 (Réf : H50C)
- Connecting kit (Réf : KRC)

### OPTIONS

- Perforated cover plate for 2-door cabinet (Réf : CAL305)
- Containment sump for 2-door cabinet (Réf : ESAL305)
- Containment sump in PVC for 2-door cabinet (Réf : SEAL305)

Model	External dimensions HxWxD (mm)	Internal dimensions HxWxD (mm)	Storage capacity (L)	Containment sump	Shelves	Load capacity	Weight (kg)
AL305	1800 x 920 x 420	1730 x 900 x 390	240		3	60	60.00