

Unicare Blue Vinyl (GS008*-A)

blue · medical grade · medium weight · vinyl examination gloves

Physical Specifications:

Material: Vinyl (latex free)
Colour: Blue
Length: 240 mm
Weight: 4.5 g
Thickness: 0.06 mm
Finish: Powder free
Fit: Ambidextrous
Sterile: No

Order Codes and Size Information:

GS0082-A | Small
GS0083-A | Medium
GS0084-A | Large
GS0085-A | Extra large

Packaging:

Each dispenser box contains 100 gloves
Each case contains 10 dispenser boxes

Classification:



AQL 1.5 (Pinholes)



EN455 (Medical Grade)



Food Safe
(not for use with fatty foods)

Features

- AQL 1.5
- Free from latex and accelerators
- Free from powder
- Tested to EN1186
- Manufactured using high-quality vinyl
- Beaded cuff
- Medium-weight construction

Benefits

- Meets medical standards for pinholes
- Suitable for those with latex allergies (Type I sensitisation) and Type IV hypersensitivity
- Reduces risks of latex allergies (Type I sensitisation) and minimises particulate contamination
- Suitable for contact with non-fatty food types
- High-quality vinyl combines low cost with reduced risk of protein allergens
- Beaded cuff provides extra strength during donning and helps prevent liquid roll back
- Medium-weight construction provides a balance of strength and sensitivity

Sectors

- Care home work
- Cleaning and janitorial
- Dentistry
- Dental hygiene
- Emergency services
- Food processing
- General assembly
- Healthcare
- Hospital and clinical duties
- Janitorial
- Local Authority work
- Painting and decorating
- Printing



For more about the Unicare range, please click on the image



UNIGLOVES UK
Registered in England No: 4010200
VAT registration 791817787

📍 3 Ambley Green, Gillingham Business Park, Gillingham, ME8 0NJ, UK
☎ +44 (0)800 049 6602
✉ enquiries@unigloves.co.uk
🌐 unigloves.co.uk

UNIGLOVES GERMANY
Arzt- & Klinikbedarf
Handelsgesellschaft mbH

📍 Camp-Spich-Straße 71
D-53842 Troisdorf-Spich,
☎ +49 (0) 2241-9323-0
✉ unigloves@unigloves.de
🌐 unigloves.de

