

MICROMIUM SERIES



STERILE

WL-CBEMS-00



Key features

- New material **MICROMIUM®**, light and very low contaminant
- Cleanroom apparel compatible from Iso class 4
- Sealed packaging by pair with an Easy Tear system
- Overpack by 25 pairs
- Apparel sterilised by gamma irradiation

Design features

- Bound seams and tunneled elastics
- Pvc Slip resistant soles
- Holding string at ankle
- Height: 50cm

Production setting

- Manufactured in a controlled environment and a facility holding ISO 9001:2008 certification

Applications

- Pharmaceutical and cosmetic industries
- Biotech and research lab
- High-tech industries
- Medical device

Packaging and tracking

100 pairs
per carton



Big bag sealed
by clip



Secondary transparent
sealed bag by 25 pairs



One pair per sealed bag
with «Easy Tear» system



Tracking number from carton to the garment

FEATURES: Overboots with PVC anti slip sole

Ties

Elastics entunneled

FABRICS AND SEAMS: Microporous film

Bound seams

PHYSICAL PERFORMANCE

Fabric Physical Properties	Test Method	Résultats	
Abrasion resistance	EN530	Class 2	
Flex cracking	ISO 7854 B	Class 6	
Trapezoidal Tear	ISO 9073-4	Class 1	
Tensile strength	ISO 13934-1	Class 1	
Puncture resistance	EN863	Class 1	

FABRIC RESISTANCE TO PENETRATION BY LIQUIDS EN ISO 6530

Fabric Chemical Properties		Penetration	Repellency
Sulphuric Acid 30%		Class 3	Class 3
Sodium Hydroxide 10%		Class 3	Class 3

STERILIZATION METHOD

GAMMA irradiation (25 kGy)

CLEANLESS

HELMKE DRUM	Test method	Particule emission rate $\geq 0.3 \mu\text{M}$	Result
	IEST-RP-CC003.3	$\leq 2\ 000$	Category 1

PACKAGING

TRIPLE PACKAGING	First	Second	Third
	2 units - Sealed and vacuum Easytear system	50 units closed by clip	200 units closed by clip

SIZE

Height	Length sole	Width sole	WeeSafe code
47 CM	36 CM	13 CM	WL-CBEMS-00

PACKING CARTON

	Units	Weight	Sizes	Volume
Carton	200	13 KG	44 X 34 X 36 CM	0,053 ^{M3}
Pallet	4 000	270 KG	80 X 120 CM	N/A

GENCOD : 3701012802459

CUSTOM CODE : 6210109800

WeeSafe – IMMEUBLE "LE FUJY"- 1 RUE DU CENTRE 93160 NOISY-LE-GRAND - FRANCE – Tél : + 33 (0)1 43 09 98 04
Mail : contact@weesafe.fr