

Product specification sheet

Revision dated 20 September 2019

1.	Product code	GM007*									
2.	Product description	Stronghold+									
3.	Material	Nitrile									
4.	Inspection parameters	Criteria		Inspection level		AQL		Standard			
		Dimensions		S-2		4.0		EN455-2			
		Physical properties		S-2		4.0		EN455-2			
		Water tight test (1,000ml)		G-1		0.65		EN455-1			
		Visual defects		G-1		2.5		In-house			
5.	Glove dimensions	Dimensions		Size (mm)							
				XS	S	M	L	XL			
		Product code	GM0071	GM0072	GM0073	GM0074	GM0075				
		Palm width	≤80	80 ± 10	95 ± 10	110 ± 10	≥ 110				
		Length	295 ± 5	295 ± 5	295 ± 5	295 ± 5	295 ± 5				
6.	Glove thickness	Point of thickness measurement			Nominal single wall (mm)						
		Cuff (at beaded end)			0.09						
		Palm (at centre of palm)			0.12						
		Fingertip (at fingertip)			0.15						
7.	Physical properties	Parameters		Before aging		After aging					
		Median force at break (N)		Min 9.0		Min 9.0					
		Tensile strength (MPa)		Min 14.0		Min 14.0					
		Ultimate elongation (%)		Min 500		Min 400					
8.	Proteins	None									
9.	Powder residue	Maximum 2.0 mg / glove									
10.	Cuff type	Beaded									
11.	Design	Ambidextrous									
12.	Surface	Microroughened									
13.	Colour	Purple									
14.	Quality assurance	Manufactured under ISO 9001 and ISO 13485									
15.	Weight (size medium)	7.3g									
16.	Standards compliant with	EN455-1	EN455-2	EN455-3	EN455-4	EN420	EN374-1	EN374-2	EN374-4	EN374-5	Food approved
		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓

The above data represents the product norm, which may be changed without notice and is intended as a guide only.