

**1. Part Number:** A316-1250

**2. Photo:**



**3. Product Information:**

Description	PharmaScoop® 1250ml
Nominal Volume:	1250ml
Material of Construction:	316L stainless steel
Method of Construction:	Crevice free, all welds ground and polished
Surface Finish:	Better than 0.5 microns Ra
Overall Length:	305mm
Body Length:	195mm
Body Diameter:	125mm
Nominal Weight:	540g
Total Surface Area:	1.24 x 10 <sup>5</sup> mm <sup>2</sup> (approx.)
BSE/TSE statement:	All polishing compounds used in the manufacture of this product are of vegetable origin
Recommended Storage Conditions:	Dry and ambient temperature

Note: This product should be thoroughly cleaned before it is used for the first time

**4. Approval:**

Signed:



Name: Kuldeep Flora  
Position: Quality Manager  
Date: 26-Oct-20