

**1. Part Number:** A316-1700

**2. Photo:**



**3. Product Information:**

|                                 |   |
|---------------------------------|---|
| Description                     | PharmaScoop® 1700ml   |
| Nominal Volume:                 | 1700ml  |
| Material of Construction:       | 316L stainless steel  |
| Method of Construction:         | Crevice free, all welds ground and polished   |
| Surface Finish:                 | Better than 0.5 microns Ra  |
| Overall Length:                 | 355mm   |
| Body Length:                    | 230mm   |
| Body Diameter:                  | 135mm   |
| Nominal Weight:                 | 700g  |
| BSE/TSE statement:              | All polishing compounds used in the manufacture of this product are of vegetable origin |
| Recommended Storage Conditions: | Dry and ambient temperature   |

Note: This product should be thoroughly cleaned before it is used for the first time

**4. Approval:**

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 26-Oct-20