


**1. Part Number:** 1210K-1500

**2. Photo:**



### 3. Product Information:

Description	Jumbo Slot Sampler 1500mm
Material of Construction:	316L stainless steel
Outer Diameter of Outer Tube:	40 mm
Inner Diameter of Inner Tube:	32 mm
Overall Length:	1550 mm
Slot Length:	120 mm
Slot Width:	20 mm
Number of Slots	4
Slot Positions:	



External Polish:	<1 Micron Ra
Internal Polish:	None
BSE/TSE statement:	All polishing compounds used in the manufacture of this product are of vegetable origin
Nominal Weight of Sampler:	4.9 Kg
Recommended Storage Conditions:	Dry and ambient temperature

### 4. Approval:

Signed:



Name: Kuldeep Flora  
Position: Quality Manager  
Date: 4-Oct-19