


1. Part Number: 8015N-170

2. Photograph:



3. Product Information:

Description:	SteriWare® Micro Spatula HDPE
Use:	Single use, disposable spatula
Material of Construction:	High Impact Polystyrene (HIPS) <ul style="list-style-type: none"> - FDA acceptable CFR 177.1640 - EU Regulations 10/2011 compliant - EC Regulations 1935/2004 compliant - BSE/TSE free
Method of Construction:	Single piece, injection moulded
Nominal Volume:	0.5ml
Overall Length:	170mm
Spatula Blade Length:	40mm
Blade Width:	6mm
Nominal Weight:	4g
Moulding & Primary Packing Environment:	Class 100,000 Cleanroom
Individually Bagged?	Yes (heat sealed PE/PET bag)
Method of Sterilisation:	N/A
Number of Scoops per Box:	100
Recommended Storage Conditions:	Dry and ambient temperature
Shelf Life:	5 years from the date of manufacture

4. Change Control:

Sampling Systems operates a policy of continual improvement.

To be kept informed of all major changes to any product in the SteriWare range please email your contact details to change.control@sampling.com

5. Country of Origin:

Made and packed in the UK.

6. Repeated Use:

This is a single-use, disposable device. Repeated use will lead to failure.

7. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 9-Mar-20