

Wilford Industrial Estate **Ruddington Lane** Wilford Nottingham **NG11 7EP**

> Tel: 0115 982 1111 Fax: 0115 982 5275

19 July 2019

TECHNICAL DATA SHEET SLE7516

Date of issue: 19.07.2019

Product Code: SLE7516

Description



30ml Universal container, polypropylene with screw cap. Product supplied in 400's, bagged in 100's.

Declaration of Conformity

We hereby confirm that this product is manufactured to the requirements of BS EN 14254:2004

QADOCPROD: SLE7516 PAGE 1 of 2 Authorised by: Andrew Littlewood





Microbiological Status	Manufactured using an automated aseptic manufacturing system under laminar flow HEPA filtration.
Standard to which the product conforms	BS EN 14254:2004
95kPa Pressure leak test to which the product conforms	To test method 49 CFR 178.605, IATA DGR 6.3.5 Ref project H20708-19 HQE2 1 Feb 2008
Material of construction (body)	Polypropylene
Cap Liner Material	Expanded low density polyethylene
Material of construction (Cap shell)	High Density polyethylene
Volume (ml)	25ml, 30ml (brim capacity), graduated to 20ml.
Bottle Dimensions (mm)	91 x 25mm (approx)
Operational temperature Range	-10 to 100 Degrees Centigrade
Inner Pack Quantity	100
Outer Box Quantity	400

QADOCPROD: SLE7516 PAGE 2 of 2, ISSUE DATE 19.07.2019

Authorised by: Andrew Littlewood



