

Material description: ACETALDEHYDE pure  
CAS No: 75-07-0

*Supplied as:* Product code Description  
CHE1000 ACETALDEHYDE pure 500ml.

## Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Purity	%	99			
Boiling range	°C	20	21		
Wt. per ml.	g			0.79	

Revision date: 09/11/1990

Issue date: 12/05/2022

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 0004/MS/19901109