

Material description: ACETIC ACID GLACIAL A.R.

CAS No: 64-19-7

<i>Supplied as:</i>	Product code	Description
	CHE1016	ACETIC ACID GLACIAL A.R. 2.5L.
	CHE1018	ACETIC ACID GLACIAL A.R. 500ml.
	CHE1020	ACETIC ACID GLACIAL A.R. 25L.

## Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Purity	%	99.8			
Boiling range	°C	117.5	119.0		
Wt. per ml.	g	1.048	1.051		
Chloride	%		0.0001		
Sulphate	%		0.0001		
Iron	%		0.00005		
Heavy metals; as Pb	%		0.00005		
Acetaldehyde	%		0.005		
Formic Acid	%		0.05		
Water	% w/w		0.2		

Revision date: 10/08/2004

Issue date: 12/05/2022

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 0011/MS/20040810