



Material description: ALUMINIUM METAL GRANULES (Coarse)

CAS No: 7429-90-5

Supplied as: Product code Description
CHE1100 ALUMINIUM METAL GRANULES (Coarse) 250g.

Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Purity	%	99.5			
Particle size	µm		600		

Revision date: 01/01/1989

Issue date: 12/05/2022

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 0033/MS/19890101