



Material description: ALUMINIUM NITRATE 9H₂O pure
CAS No: 7784-27-2

Supplied as: Product code Description
CHE1104 ALUMINIUM NITRATE 9H₂O pure 250g.

Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Purity	%	98			

Revision date: 01/01/1989

Issue date: 12/05/2022

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 0795/MS/19890101