

Material description: AMMONIUM FERROUS SULPHATE 6H₂O pure
CAS No: 7783-85-9

<i>Supplied as:</i>	Product code	Description
	CHE1166	AMMONIUM FERROUS SULPHATE 6H ₂ O pure 500g.
	CHE1168	AMMONIUM FERROUS SULPHATE 6H ₂ O pure 250g.
	CHE1170	AMMONIUM FERROUS SULPHATE 6H ₂ O pure 1Kg.

Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Purity	%	98			
Copper	%		0.005		
Zinc	%		0.02		

Revision date: 01/01/1989

Issue date: 12/05/2022

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 0834/MS/19890101