



Material description: AMMONIUM OXALATE H2O pure
CAS No: 6009-70-7

<i>Supplied as:</i>	Product code	Description
	CHE1196	AMMONIUM OXALATE H2O pure 100g.
	CHE1198	AMMONIUM OXALATE H2O pure 250g.

Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Purity	%	98.0	101.0		

Revision date: 22/11/1999

Issue date: 12/05/2022

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 0071/MS/19991122