

Material description: CYCLOHEXANONE pure
CAS No: 108-94-1

| | | |
|---------------------|--------------|---------------------------|
| <i>Supplied as:</i> | Product code | Description |
| | CHE1742 | CYCLOHEXANONE pure 500ml. |
| | CHE1744 | CYCLOHEXANONE pure 2.5L. |

Specification

| Parameter | Units | Minimum value | Maximum value | Typical value | Notes |
|---------------|-------|---------------|---------------|---------------|-------|
| Purity | % | 99 | | | |
| Boiling range | °C | 153 | 156 | | |
| Wt. per ml. | g | | | 0.94 | |

Revision date: 24/09/1991

Issue date: 12/05/2022

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 1006/MS/19910924