

Material description: DICHLOROMETHANE pure  
CAS No: 75-09-2

<i>Supplied as:</i>	Product code	Description
	CHE1828	DICHLOROMETHANE pure 2.5L.
	CHE1830	DICHLOROMETHANE pure 500ml.
	CHE1832	DICHLOROMETHANE pure 25L.

## Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Boiling range	°C	39.0	40.5		
Wt. per ml.	g	1.323	1.327		
Water	% w/w		0.05		
Non-volatiles	%		0.01		

Revision date: 09/05/2018

Issue date: 12/05/2022

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 0237/MS/20180509