

Material description: DIETHYL ETHER A.R.

CAS No: 60-29-7

<i>Supplied as:</i>	Product code	Description
	CHE1836	DIETHYL ETHER A.R. 2.5L.
	CHE1838	DIETHYL ETHER A.R. 500ml.

## Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Purity	%	99.5			
Boiling range	°C	34	36		
Wt. per ml.	g	0.712	0.716		
Water	% w/v		0.02		
Non-volatiles	%		0.002		

Revision date: 03/07/2020

Issue date: 12/05/2022

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 0251/MS/20200703