

Material description: ETHYL BENZOATE pure
CAS No: 93-89-0

Supplied as: Product code Description
CHE1970 ETHYL BENZOATE pure 100ml.

Specification

| Parameter | Units | Minimum value | Maximum value | Typical value | Notes |
|---------------|-------|---------------|---------------|---------------|-------|
| Purity | % | 99 | | | |
| Boiling range | °C | 211 | 213 | | |
| Wt. per ml. | g | 1.051 | 1.052 | | |

Revision date: 01/01/1989 Issue date: 12/05/2022

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 1097/MS/19890101