

Material description: FEHLING'S SOLUTION NO. 2 QUALITATIVE
CAS No: Mixture

Supplied as: Product code Description
 CHE1984 FEHLING'S SOLUTION NO. 2 QUALITATIVE 1L.
 CHE1986 FEHLING'S SOLUTION NO. 2 QUALITATIVE 500ml.

Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Functional test <i>by Fehling's test</i>					Red precipitate formed
Density @20°C	g/ml	1.248	1.268	1.258	

Revision date: 30/04/2020 Issue date: 12/05/2022

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 1120/MS/20200430