

Material description: IRON(II) SULPHATE 7H<sub>2</sub>O A.R.

CAS No: 7782-63-0

Supplied as: Product code Description  
CHE2014 IRON(II) SULPHATE 7H<sub>2</sub>O A.R. 500g.

## Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Purity	%	99.5			
Lead	%		0.0002		
Arsenic	%		0.0001		
Cadmium	%		0.0001		
Mercury	%		0.00001		

Revision date: 19/05/2015

Issue date: 12/05/2022

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 1535/MS/20150519