

Material description: HYDROCHLORIC ACID 2.0M (2N)

CAS No: 7647-01-0

<i>Supplied as:</i>	Product code	Description
	CHE2124	HYDROCHLORIC ACID 2.0M (2N) 1L.
	CHE2126	HYDROCHLORIC ACID 2.0M (2N) 2.5L.
	CHE2128	HYDROCHLORIC ACID 2.0M (2N) 5L.

Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Strength	mol/L	1.998	2.002		

Revision date: 01/01/1989

Issue date: 12/05/2022

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 1734/MS/19890101