



**SCIENTIFIC  
LABORATORY  
SUPPLIES**

Material description: n-BUTYRALDEHYDE pure  
CAS No: 123-72-8

Supplied as: Product code Description  
CHE2616 n-BUTYRALDEHYDE pure 250ml.

## Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Purity	%	98			
Boiling range	°C	74	78		
Wt. per ml.	g	0.802	0.804		

Revision date: 17/05/1990 Issue date: 18/09/2025

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 1685/MS/19900517