



**SCIENTIFIC
LABORATORY
SUPPLIES**

Material description: PETROLEUM SPIRIT 60/80 pure

CAS No: 8032-32-4

<i>Supplied as:</i>	Product code	Description
	CHE2806	PETROLEUM SPIRIT 60/80 pure 500ml.
	CHE2808	PETROLEUM SPIRIT 60/80 pure 1L.
	CHE2810	PETROLEUM SPIRIT 60/80 pure 2.5L.

Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Boiling range; 90%	°C	60	80		
Wt. per ml.	g	0.660	0.680		

Revision date: 11/07/1990

Issue date: 12/05/2022

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 1303/MS/19900711