



**SCIENTIFIC
LABORATORY
SUPPLIES**

Material description: PHENOLPHTHALEIN SOLUTION

CAS No: Mixture

<i>Supplied as:</i>	Product code	Description
	CHE2846	PHENOLPHTHALEIN SOLUTION 500ml.
	CHE2848	PHENOLPHTHALEIN SOLUTION 100ml.
	CHE2850	PHENOLPHTHALEIN SOLUTION 2.5L.

Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Colour change; Clear colourless to Red					pH 8.3-10.0

Revision date: 01/01/1989

Issue date: 12/05/2022

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 1308/MS/19890101